

Dermat+

| Patch Testing Catalogue



CONTENTS

Hapten series

Patch test chambers
& Skinmarkers

Spot tests

Hapten series

Expiry & stock keeping	1
1.1 Belgian Standard Serie	
1.1.1 Belgian core baseline serie	2 - 3
1.1.2 Belgian extended baseline serie	3
1.1.3 Belgian supplements (BCEDG)	4
1.2 Belgian Cosmetic & Pharmaceutical Series	
2.1.1 Belgian Cosmetic Serie	5
2.1.2 Belgian Pharmaceutical Serie	5
1.3 International Screening Series	6

Expiry & stock keeping

Hapten

Storage

Store at 2° - 8°C*

During transport haptens are subject to ambient temperatures. Exposure to ambient temperatures for haptens during transit (< 10 days) will not affect product stability* as this is insufficient time for influence homogeneity or concentration. Long term storage of haptens should however always be refrigerated.

Expiry

30, 24, 12 or 8 months

Haptens are divided into one out of five different expiry groups depending on their composition. As all haptens must be delivered to the end user with a pre defined time remaining until expiration, planning of order volumes and logistics must take expiration dates into account. The expiry groups are defined by how many months the hapten is stable from the date of production.

The expiry groups are: 30, 24, 18, 12 and 8 months.

30 and 24 months

The 30 and 24 months haptens must all be delivered to the end user with a minimum of 12 months remaining until expiration. Haptens in these groups are guaranteed to have a shelf life of at least 14 months when they are dispatched from the Chemotechnique premises.

12 months

The 12 months haptens must all be delivered to the end user with a minimum of 6 months remaining until expiration. Haptens in these groups are guaranteed to have a shelf life of at least 8 months when they are dispatched from the Chemotechnique premises.

8 months

The 8 months haptens must all be delivered to the end-user with a minimum of 6 months remaining until expiration. As Chemotechnique may stock these haptens for up to two months before dispatch, the 8 month haptens may never be kept in stock at the distributor and should therefore only be ordered when requested by end users.

Hapten series



1.1 BELGIAN STANDARD SERIE



Modifications

1.1.1 BELGIAN CORE BASELINE SERIE

		* expiry (months)	
1.	Potassium dichromate 0.5% pet	*30 m	P-014A
2.	Neomycin sulfate 20% pet	*30 m	N-001
3.	Thiuram mix 1% pet	*30 m	Mx-01
4.	P-phenylenediamine (PPD) 1% pet	*30 m	P-006
5.	Cobalt chloride hexahydrate 1% pet	*30 m	C-017A
6.	** Caine III Mix 10% pet	*30 m	Mx-19
7.	Formaldehyde 2% aqua	*12 m	F-002B
8.	Colophonium 20% pet	*30 m	C-020
9.	Propolis 10% pet	*30 m	P-022
10.	Balsam of Peru 25% pet	*30 m	B-001
11.	N-Isopropyl-N-phenyl-paraphenylenediamine (IPPD) 0,1% pet	*30 m	I-004
12.	Wool alcohols = Lanolin alcohol 30% pet	*30 m	W-001
13.	Mercapto mix 2% pet	*30 m	Mx-05A
14.	Epoxy resin, bisphenol A 1% pet	*30 m	E-002
15.	Paraben mix 16% pet	*30 m	Mx-03C
16.	4-tert-Butylphenolformaldehyde resin (PTBP) 1% pet	*30 m	B-024
17.	Fragrance mix I 8% pet	*30 m	Mx-07
18.	Sesquiterpene Lactone mix 0,1% pet	*30 m	Mx-18
19.	**2-Hydroxyethyl methacrylate 2% pet	*12 m	H-010
20.	Nickel(II)sulphate hexahydrate 5% pet	*30 m	N-002A
21.	*Methylisothiazolinone+Methylchlorisothiazolinone 0,02% aqua	*12 m	C-009B
22.	Mercaptobenzothiazole (MBT) 2% pet	*30 m	M-003A
23.	Tixocortol-21-pivalate 0.1% pet	*30 m	T-031B



* expiry (months)

24.	Budesonide 0,01% pet	*30 m	B-033B
25.	Methylidibromo Glutaronitrile 0,5% pet	*30 m	D-049E
26.	Hydroxyisohexyl 3-cyclohexene carboxaldehyde (Lyrall) 5% pet	*12 m	L-003
27.	Fragrance Mix II 14,0% pet	*30 m	Mx-25
28.	*Methylisothiazolinone 0.2% aqua	*12 m	M-035B
29.	Textile Dye Mix 6.6% pet	*30 m	DC342
30.	Evernia Furfuracea 1% pet	*30 m	E-026
31.	Sodium Metabisulfiet 1% pet	*30 m	S-011
32.	Benzisothiazolinone 0.1% pet	*30 m	B-003B
33.	Decyl glucoside 5% pet	*30 m	D-065

1.1.2 BELGIAN EXTENDED BASELINE SERIE

1.	2-BROMO-2-NITROPROPANE-1,3-DIOL 0,5% pet	*24 m	B-015B
2.	Diazolidinyl Urea (Germall II) 2% pet	*30 m	D-044A
3.	Compositae mix II 5 % pet	*30 m	Mx-29A
4.	Parthenolide 0.1% pet	*30 m	P-029
5.	Hydroperoxides of Limonene 0.3% pet	*8 m	H-032A
6.	Hydroperoxides of Limonene 0.2% pet	*8 m	H-032B
7.	Hydroperoxides of Linalool 1% pet	*8 m	H-031A
8.	Hydroperoxides of Linalool 0.5% pet	*8 m	H-031B
9.	2-n-Octyl-4-isothiazolin-3-one 0.1% pet	*30 m	O-004
10.	Sorbitan sesquioleate 20% pet	*30 m	S-005
11.	Sorbitan oleate 5% pet	*30 m	S-004



* expiry (months)

1.1.3 BELGIAN SUPPLEMENTS (BCEDG)

1.	Carvone 5% pet	*30 m	C-035
2.	Propylene glycol 30% aqua	*12 m	P-019B
3.	Diphenylguanidine 1% pet	*30 m	D-022
4.	Black Rubber Mix 0.6% pet	*30 m	Mx-04
5.	Methyldibromo glutaronitrile 0,3% pet	*30 cm	D-049A
6.	Amerchol L - 10150% pet	*30 cm	A-004
7.	Aluminium (III)chloride hexahydrate 2% pet	*30 cm	A-022

OPGELET

- ** Soms zwakke erythemateuze irriterende reacties : Een 2de lezing na 7 dagen is aangewezen.
- * Gebruik een micro-pipet! 15 µl voor Finn-kamers, 20 µl voor IQ Ultra-kamers; gebruik anders lagere concentraties om mogelijke overgevoeligheid te voorkomen. Vals-negatieve reacties zijn niet uit te sluiten.

ATTENTION

- ** Parfois de légères réactions d'irritations érythémateuses : Une deuxième lecture après 7 jours est recommandée.
- * Utiliser une micro-pipette ! 15 µl pour les chambres Finn, 20 µl pour les chambres IQ Ultra ; sinon, utiliser des concentrations plus faibles pour éviter une hypersensibilité. Toutefois des réactions fausses/ négatives peuvent apparaître.

! I-001A / C-007 / L-004 / L-002A behoren niet langer tot de Belgische Standaard reeks.

! I-001A / C-007 / L-004 / L-002A ne font plus partie de la batterie standard belge.



1.2 BELGIAN COSMETIC AND PHARMACEUTICAL SERIES

1.2.1 BELGIAN COSMETIC SERIE

* expiry (months)

1.	Cetrimide 0.05% aqua	*30 m	AP996
2.	Benzoic Acid 5% pet	*30 m	B-005
3.	Sorbic Acid 2% pet	*30 m	S-003
4.	Dexpanthenol 5% pet	*30 m	P-042 / MS373
5.	Cetearyl Alcohol 20% pet	*30 m	C-033
6.	Phenoxyethanol 1% pet	*30 m	P-025
7.	Cocamidopropylbetaine 1% aqua	*12 m	C-018
8.	Dimethylaminopropylamine 1% aqua	*12 m	D-053
9.	Iodopropynyl butylcarbamate 0.2% pet	*30 m	I-008C
10.	Propyl gallate 1% pet	*24 m	P-021
11.	Ethylhexylglycerin 5% pet	*24 m	E-027
12.	Chlorhexidine digluconate 0.5% aqua	*12 m	C-005

1.2.2 BELGIAN PHARMACEUTICAL SERIE

1.	Chloramphenicol 5% pet	*30 m	C-032
2.	Nitrofurazone 1% pet	*30 m	N-005
3.	Fusidic Acid (Sodium Salt) 2% pet	*30 m	F-003
4.	P-Chloro-m-Cresol (PCMC) 1% pet	*24 m	C-008
5.	Benzoyl peroxide 1% pet	*12 m	B-007
6.	Isobetadine (pure)		

Gelieve ons te contacteren voor alle bestellingen van International Screening Series.
Veuillez nous contacter pour toutes commandes de International Screening Series.

1.3 INTERNATIONAL SCREENING SERIES

1.	Bakery	B-1000
2.	Corticosteroid	CS-1000
3.	Cosmetic	C-1000
4.	Cutaneous Adverse Drug Reaction	CAD-1000
5.	Dental Screening	DS-1000
6.	Dental Materials - Patients	DMP-1000
7.	Dental Materials - Staff	DMS-1000
8.	Epoxy	E-1000
9.	Fragrance	F-1000
10.	Hairdressing	H-1000
11.	Implant	IMP-1000
12.	Isocyanate	I-1000
13.	Leg Ulcer	LU-1000
14.	Medicament	ME-1000
15.	Metal	MET-1000
16.	Metal Extended	METE-1000
17.	(Meth) Acrylate - Adhesives, Dental, Printing & Other	MA-1000
18.	(Meth) Acrylate - Artificial Nails	MN-1000
19.	Oil & Cooling Fluid	O-1000
20.	Photopatch	PP-1000
21.	Plant	PL-1000
22.	Plastic & Glue	PG-1000
23.	Rubber Additives	R-1000
24.	Shoe	SH-1000
25.	Sunscreen	SU - 1000
26.	Textile Colors & Finish	TF-1000

Patch test chambers & Skinmarkers

2.1 Patch test chambers

2.1.1 Patches	8
2.1.2 Transparant films	9
2.1.3 Accessories	9

2.2 Skinmarkers

2.2.1 Sterile	10
2.2.2 Non-sterile	10

Patch test chambers & Skinmarkers



2.1 PATCH TEST CHAMBERS

2.1.1 PATCHES



BasiQ ULTRA - Patch test unit

BIQ-U	Comfortable	1 box
	Hypoallergenic - latex free	50 x 10 pc.



IQ ULTRA

IQ-U	Comfortable	1 box
	Hypoallergenic - latex free	100 x 10 pc.



ALLERGEAZE

4300	Comfortable	1 box
	Hypoallergenic - latex free	100 x 10 pc.



VAN DER BEND

4600	Comfortable	1 box
	Hypoallergenic - latex free	50 x 10 pc.
	non woven	



FINN CHAMBERS

10000010	Round aluminium chambers	1 box
	hypoallergenic - latex free	100 x 10 pc.



CURATEST (L&R)

163055	Epicutaneous test patch	1 box
	of non woven fabric	50 x 10 pc.
	hypoallergenic - latex free	



2.1.2 TRANSPARENT FILMS



IQ ULTIMATE


IQ-UL	Elastic, transparent, water-resistant Hypoallergenic - latex-free	1 box 100 x 10 pcs
-------	---	-----------------------



ALLERGEAZE CLEAR

ALLERGEAZE -CLEAR	Elastic, transparent, water-resistant Hypoallergenic - latex-free	1 box 100 x 10 pcs
----------------------	---	-----------------------

2.1.3 ACCESSORIES



Application device (43 x 27cm)

AP - P	Helpful tool to prepare the ULTIMATE and IQ ULTRA patch chambers.
--------	--



Reading plate (IQ ULTRA & IQ ULTIMATE)

RP -P	Facilitates the reading of the patch test results.
-------	--



Lid - tray for chambers - Allergeaze

23SPC-1003



Tray for chambers - Allergeaze

23SPC-1004



2.2 SKINMARKERS

2.1.1 NON - STERILE



Chemo Skinmarker Medium

SMM-4 Gentian violet ink
Round tip
4 pcs / box

2.1.2 STERILE



Skinmarker Sterile

19005-02CE Gentian violet ink
Fine tip
Available with flexible ruler

Spot tests

3.1 Chemo Nickel Spot test	12
3.2 Chemo Cobalt Spot test	12

Spot tests



3.1 CHEMO NICKEL SPOT TEST



NT

A simple test to detect nickel in metal objects. The test detects free nickel down to a limit of 10 ppm (parts/million).

3.2 CHEMO COBALT SPOT TEST



CT CHEMO

A simple test to detect cobalt in metal objects. The test detects free cobalt down to a limit of 8.3 ppm (parts/million).

Wenst u meer info over reeksen en testen?

Vraag naar de interactieve pdf

via info@dermat.be

Vous souhaitez plus d'infos sur les batteries et les tests?

Demandez le pdf interactif à info@dermat.be

Dermat Medical Supplies

Ambachtenlaan 38

3001 Leuven

info@dermat.be

016 40 05 45

www.dermat.be



Dermat+