

medisun® PDT 1200



Derma^t

info@dermat.be
www.dermat.be

THE IDEAL ADDITION FOR THE EXPOSURE OF
SMALL-TO MEDIUM-SIZED AREAS ALL OVER
THE BODY

ADVANTAGES MEDISUN® PDT 1200

- Compact and mobile PDT- System
- Microcontroller control with dose input in J/cm²
- With Multi-Colour-PC-Control (as option), integrated 3D- fluorescence diagnosis (woodlight), red-light, blue-light, SAD-therapy, individually programmable light spectra and therapy protocols
- Easy positioning over 3 axes
- Flexible and space saving
- LED lamp field with 12 high-performance-LEDs
- Individual adjustment of parameters – future-proof and flexible
- The exposure is painless. During the exposure the patient can listen to music, read or just relax
- Only 1 hour exposure time needed with the clinical approved painless standard protocol
- Therapy is year-round predictable and can always be performed in a controlled way with a homogeneous and reproducible dose
- Small- to medium-sized lesions can be treated individually on the whole body
- The therapy is independent of the weather
- Only a regular 230 V-plug is needed
- Low operating costs due to the latest LED technology
- Comparatively low investment costs will be amortized very quickly, usually within a few months
- You can achieve additional income and gain many satisfied patients
- All medisun® devices are developed and manufactured in Germany
- Made in Germany – CE-certified according to MDD 93/42/EEC



THERAPEUTIC OPTIONS:

- Actinic Keratosis
- Skin Rejuvenation
- And other indications....



SCHULZE & BÖHM GmbH
Hermülheimer Str. 10
D-50321 Brühl
Germany

Tel.: +49 (0) 2232-57998-0
Fax: +49 (0) 2232-57998-11

E-Mail: info@medisun.de
www.medisun.de